



# POWER SOLUTIONS

## for Patient Care Areas



### PowerMATE® Available 15 amp or 20 amp

- 15 amp circuit breaker or 20 amp breaker
- Mounts to an IV pole with a diameter of 0.875" to 1.50"
- 15 foot hospital grade power cord
- UL recognized component to UL1363A
- 1-year warranty
- 15 amp - PS11076 | 20 amp - PS12595

### PowerMATE® - TR

- Tamper Resistant Safety Outlets
- 15 amp breaker
- Mounts to an IV pole with a diameter of 0.875" to 1.50"
- 15 foot hospital grade power cord
- UL recognized component to UL1363A
- 1-year warranty
- PS12744

**NEW!**



### IsoMATE®

- 5 amp breaker
- Fully contained isolation transformer offers full line isolation
- Less than 100 microamps leakage current
- Dual circuit breakers provide overload protection
- Mounts to an IV pole with a diameter of 0.75" to 2.50"
- Regulatory Approvals: UL listed to 60601-1, CSA C222 No. 601
- 2-year warranty
- PS11736

Call us at 888.702.9686 or email [aivsales@aiv-inc.com](mailto:aivsales@aiv-inc.com)  
[www.aiv-inc.com](http://www.aiv-inc.com)